

NMR

5mm Sample tubes

- ▶ Ultra Precision
- ▶ Factory Pricing
- ▶ Mass Production



NES-600 for ~600MHz

- Length: 7inch(178mm)
- Wall thickness:0.38mm

NEW!
Quality from JAPAN
Request
a Free Sample!



• 1H spectrum with sample spinning (20Hz), NS=16
(Sample : Methyl-d3 Alcohol-d)



Model	Outer Dia.	Concentricity	Camber	Range	Package
NGS-400	4.950+0.038 mm	±0.01 mm	±0.06 mm	~ 400 MHz	50ea / pack
NES-600	4.950±0.018 mm	±0.01 mm	±0.025mm	~ 600 MHz	10ea / pack

www.optimajp.com

A Technomart Company / An ISO 9001 Certified / tmk@technomartinc.com