

NMR



5 mm Sample tubes

- ▶ Ultra Precision
- ▶ Factory Pricing
- ▶ Mass Production

NLS-800 for ~800MHz

- Length: 7inch(178mm)
- Wall thickness:0.38mm

NEW!
Quality from JAPAN
Request
a Free Sample!



Model	Outer Dia.	Concentricity	Camber	Range	Package
NLS-800	4.956 ± 0.008 mm	± 0.004 mm	± 0.015 mm	~800MHz	10ea / pack
NTS-1000	4.956 ± 0.008 mm	± 0.0025 mm	± 0.0075 mm	~950MHz	10ea / pack

MADE IN JAPAN

www.optimajp.com

A Technomart Company / An ISO 9001 Certified / trmk@technomartinc.com