

MICRO CENTRIFUGE

- **Compact and Noise Free Design**

Small foot print conserves valuable bench space.

- **Max. 10,000RPM**

Designed to operate at 10,000rpm. Ideal for quick spin-downs of sample tube walls or cap.

- **8-place Rotor**

Fixed-angle rotor for eight 1.5ml micro centrifuge tubes, or 0.4ml tubes with adaptors provided.

- **Neon Lighting**

Unique see-through design with attractive neon lighting.

- **Stability**

Specially designed rubber feet prevent the unit from moving during operation.



OS-100B Micro Centrifuge

Specifications

Model	OS-100B
Speed	10,000RPM
G-force	5,590 x g
Capacity	1.5ml x 8 tubes
Drive	Direct drive system
Dimensions	153 x 162 x 122 (H)mm
Lighting	Neon
Weight	1kg
Electrical	220V, 60Hz / 100 - 110V, 50/60Hz



Optional Accessory

Model	Description
V100B1	Adapter for 0.4ml, 8ea/set



- * Continuing research and improvements may result in specification change at any time.
- * 1 year warranty on all parts.
- * CE marked. Copies of CE documentation are available for compliance verification.



World Headquarters

OPTIMA INC.

Rogie Bldg. 5F, 1-48-11 Itabashi, Itabashi-ku,
Tokyo 173-0004 Japan
Tel.+81-3-5375-2351 Fax.+81-3-5375-2360
E-mail: optima@optimajp.com

U.S.A Operations

OPTIMA USA, INC.

2254-B Landmeier Road, Elk Grove Village
IL 60007 U.S.A
Tel.+1-847-758-8834 Fax.+1-847-758-8836
E-mail: optima@optimajp.com

Distributor